

# WEB322 - COMT expression

Benjamin Weigel

17th November, 2014

## 1 Day 1

- 500 mL over night culture (+ 200  $\frac{\mu g}{ml}$  Amp) of **ZYP-5052** (autoinduction medium) inoculated with toothpick of COMT (3) from Th. Vogt

## 2 Day 2

- measured  $OD^{600} = 10.56 \rightarrow 500 \mu l$  **SDS-PAGE** sample
- harvested cells by centrifugation (10.000 x g, 4°C , 5 min)
- inoculated 250 mL culture of **LB** (+ 200  $\frac{\mu g}{ml}$  Amp) with 1.8 mL of overnight AI culture
- started incubation at 10.50

time	sample	$OD^{600}$	comments
10.50	AI	10.56	
13.00	LB culture	0.852	
18.00	LB culture	3.6	SDS-PAGE sample NI

## 3 Purification from AI-medium expression

- resuspended pellet ( $\sim 3$  g) in 60 mL **Lysis Buffer** (50 mM Tris/HCl, 0.5 M NaCl, 10% (v/v) glycerol, 10 mM Imidazole, 1% Tween 20)
- added a spatula tip of HEWL and incubated for 15 min at RT while inverting
- lysed cells by sonication (3 times for 30 s at 70% amplitude, 1s-1s on-off-cycle)
- added DNase buffer and 200  $\mu l$  DNase I and incubated on ice for 15 min
- removed debris by centrifugation (20 min, 10.000 x g, 4°C )

### 3.1 ÄKTA

- filtered through a 0.45  $\mu\text{m}$  filter and injected onto 1 mL HiTrap Talon FF column on ÄKTA
- ! buffers, lysate and also the collection tubes were cooled
- column was equilibrated with **buffer A** (50 mM Tris/HCl, 500 mM NaCl, 10% (v/v) glycerol, 2.5 mM imidazole, pH 7.4)
- after sample injection the column was washed with 4 CV of buffer A and 4 CV of 5% buffer B
- eluted with 8 CV **buffer B** (50 mM Tris/HCl, 500 mM NaCl, 10% (v/v) glycerol, 300 mM imidazole, pH 7.4)
- !!! only a very small protein peak elutes at 100% B (Fig. 1a)  $\rightarrow$  SDS-Gel to clarify what went wrong in purification (Fig. 2a)
- precipitated 100  $\mu\text{l}$  of fractions X1 and A3 with TCA and used for SDS-PAGE

## 4 Purification from LB-medium expression

- harvested cells by centrifugation (10.000 x g, 4°C , 5 min) and stored at -20°C






Figure 2: COMT purification.

**2a:** COMT purification using 1 mL HiTrap Talon FF (left side). NI - after induction (total), S - after induction (soluble), X - not sure (discard), FT - ÄKTA flowthrough, X1 - fraction X1, A3 - elution fraction A3.

**2b:** Flowthrough, that was collected during COMT purification using 1 mL HiTrap Talon FF. The flowthrough stood at room temp for ~3 hours. The turbid solution was centrifuged to separate soluble and aggregated protein. S - soluble, I - insoluble.

**2c:** Flowthrough, that was collected during first COMT purification trial (2a) was diluted with 3 volumes of **buffer A.2** (50 mM Tris/HCl, 500 mM NaCl, 10% (v/v) glycerol, pH 7.4). 30 mL of the diluted flowthrough (2.5 mM imidazol concentration) were applied to a HisTrap (Ni-NTA) 1 mL FF column at 0.5 mL/min flow. X1 - flowthrough, A1:4 - fractions A1 to A4.