

Protocol for electrocompetent Agro cells

Glycerol stocks:
for GV3101 (main strain; +Rifampicin)

ON culture + Rifampicin (28°C)

Add ON culture in the main culture (+ Rifampicin) and incubate for 15 – 18 hours (overnight) at 28°C

Divide culture in big centrifuge tubes

All steps on ice

Centrifuge: 3000xg, 10 min, 4°C

Wash 2 times with 10% Glycerol (cool down before using)

After the second wash step and third centrifugation mix the cells together

Centrifuge and resuspend the pellet in 1/800 10% Glycerol

250 µL Glycerol for 200 mL culture

625 µL Glycerol for 500 mL culture

1,25 mL Glycerol for 1 L culture

1,5 mL Glycerol for 1,2 L

aliquot 25 µL/tube, shock frost in liquid nitrogen

store at -80°C

ctal

107

A

B

1.72

pBEW 107

232.5

9.768 ml

1.32

control

303

1.697 ml

Glycerol		autoclave always
10%	- 126g	1L
	- 252g	2L
	- 378g	3L
30%	- 378g	1L
	- 189g	500mL
	- 94,90g	250mL

Trafo electrocomp. Agros
 (1 μ L) DNA 0,5 μ L also ok
 + 25 μ L cells GV3101 (on ice)
 add to electroporation cuvette
 → pulse
 add 1mL LB
 pipette to tube and shake
 at 28°C 3-4 hours
 centrifuge and resuspend for plating

Dissolving buffer / Infiltration

100 mL	100 mM MES (pH 5.5)
10 mL	1 M MgSO_4
<u>1 L</u>	

100 mM MES - 1 L - 19,52 g pH 5.5 NaOH
 1 M MgSO_4 - 100 mL - 24,65 g

10^{-1} Verdünnung AgtOS

$\frac{0,2}{\text{OD Wert } 600}$ • 2000 \rightarrow $\pm \mu\text{L}$ + Rest
 \rightarrow gesamt vol. (total 2 mL)

\rightarrow o/n cult. 1:10 mit H_2O verd.
 für die Rechnung Wert
 wieder $\times 10$ nehmen

10^{-2} Verdünnung AgtOS

$\frac{0,02}{\text{OD Wert } 600}$ • 2000 \rightarrow $\pm \mu\text{L}$ + Rest
 \rightarrow gesamt vol. (total 2 mL)

\rightarrow o/n cult. 1:10 mit H_2O verd.
 für die Rechnung aber Wert
 wieder $\times 10$ nehmen

Verdünnungen f. Infiltrationen

1:10 $\rightarrow 10^{-1}$	200 μL / 2 mL
1:100 $\rightarrow 10^{-2}$	20 μL / 2 mL
1:1000 $\rightarrow 10^{-3}$	2 μL / 2 mL
1:10.000 $\rightarrow 10^{-4}$	0,2 μL / 2 mL
1:100.000 $\rightarrow 10^{-5}$	0,02 μL / 2 mL

oder aus
 Verdünnungen \rightarrow

Bsp. 10^{-2} Verd. 20 μL / 2 mL

\sim 200 μL 10^{-2} Verd. / 2 mL $\rightarrow 10^{-3}$

\sim 200 μL 10^{-3} Verd. / 2 mL $\rightarrow 10^{-4}$

\sim 200 μL 10^{-4} Verd. / 2 mL $\rightarrow 10^{-5}$

\sim 200 μL 10^{-5} Verd. / 2 mL $\rightarrow 10^{-6}$
 usw.