

WEB305 - Chromabond HR-X Test

Benjamin Weigel

August 27, 2014

1 Introduction

For the analysis of methyl transferase reaction via the LAAO-assay it is necessary to remove flavonoids and phenyl propanoids. Those substances interfere with the detection of H_2O_2 via HRP, because rather than TMB these substances tend to get oxidized.

Trying to remove the flavonoid substrates and products of a MT-reaction by solid phase extraction columns(SPE).




Figure 1: Substances to be separated for analysis by LAAO assay. Polar substances to be eluted for analysis are **1** - SAM, **2** - SAH and **3** - Ribosyl homocysteine. Flavonoids and phenyl propanoids (**4–5**) are problematic during detection with LAAO/HRP.

2 Experimental

2.1 Sample preparation

Sample solutions:

1. 100 μM SAH in ddH₂O
2. 100 μM SAH, 0.2 mM Eriodictyol in ddH₂O

Sample treatment:

- A1:** 100 μM SAH in ddH₂O passed through HR-X column
A2: 100 μM SAH, 0.2 mM Eriodictyol in ddH₂O passed through HR-X column
B1: 100 μM SAH in ddH₂O
B2: 100 μM SAH, 0.2 mM Eriodictyol in ddH₂O

Column preparation:

Conditioning of HR-X column:

1.5 mL MeOH, then 1.5 mL Water

Sample aspiration:

0.5 mL sample passed through column by pressure (syringe)

Washing/Elution:

eluted polar components by washing with 0.5 mL water

Regeneration:

regenerated by passing through 2 mL of MeOH

2.2 LAAO-Reaction

1. add 75 μ L LAAO-reagent to 25 μ L of sample
2. incubate for 10 min @ RT
3. stop reaction with 50 μ L 12 N H₂SO₄
4. measure absorption at 450 nm