

# WEB326 - COMT expression

Benjamin Weigel

25th November, 2014

## 1 Day 1

- 500 mL over night culture (+ 200  $\frac{\mu g}{ml}$  Amp) of **ZYP-5052** (autoinduction medium) inoculated with toothpick of COMT (culture that stood at 4°C for 2 weeks)

## 2 Day 2

- measured OD<sup>600</sup> =  $\sim$  2.6 at 10.00
- measured OD<sup>600</sup> = 6.88  $\rightarrow$  at 1300
- harvested cells by centrifugation (10.000 x g, 4°C , 3 min)

!! cells grew slowly  $\rightarrow$  probably due to old pre-culture....

## 3 Purification from AI-medium expression

- resuspended pellet ( $\sim$ 3 g) in 60 mL **Lysis Buffer** (50 mM Tris/HCl, 0.5 M NaCl, 10% (v/v) glycerol, 10 mM Imidazole, 1% Tween 20)
- added a spatula tip of HEWL and incubated for 15 min at RT while inverting
- lysed cells by sonication (3 times for 30 s at 70% amplitude, 1s-1s on-off-cycle)
- added DNase buffer and 200  $\mu$ l DNase I and incubated on ice for 15 min
- removed debris by centrifugation (20 min, 10.000 x g, 4°C )

### 3.1 ÄKTA

- filtered through a 0.45  $\mu$ m filter and injected onto 1 mL HiTrap Talon FF column on ÄKTA
- ! buffers, lysate and also the collection tubes were cooled
- column was equilibrated with **buffer A** (50 mM Tris/HCl, 500 mM NaCl, 10% (v/v) glycerol, pH 7.4)
- after sample injection the column was washed with 8 CV of buffer A, 8 CV of 5% buffer B and 8 CV of 10% buffer B

- eluted with 8 CV **buffer B** (50 mM Tris/HCl, 500 mM NaCl, 10% (v/v) glycerol, 300 mM imidazole, pH 7.4)

!!! only a very small protein peak elutes at 100% B (Fig. ??) → SDS-Gel to clarify what went wrong in purification