
24.09.2012 11:12:41

Enhancement: 1 Brightness: 65
Exposure Time: 1774 Contrast: 52
Gamma: 36

Colony PCR
pET28a (P)ERT + phoB
Doppeltransformante

O:\ABT\NWC\Weigel_WEB\BioDoc\120828_Luxs&SA
HH.TIF

06.12.2012 09:37:47


Enhancement: 1
Exposure Time: 45

Brightness: 50
Contrast: 50
Gamma: 25

Gamma: 25

Muscle - stiffness
Contractant
Contractile
Contractile tissue
Post-tend.

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28.08.2012 16:30:40

Enhancement: 27 Brightness: 50
Exposure Time: 50 Contrast: 50
Gamma: 25

27.07.2012 14:24:25

Enhancement: 5 Brightness: 57
Exposure Time: 40 Contrast: 74
Gamma: 25

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
Fig. 1. The nucleotide sequence of *E. coli* MTA /SAH nucleosidase annotated with the deduced amino acids. The open reading frame was amplified as an *EcoRI/NotI* fragment using the primer set: 5'-CTC **CAA** **TTC** TCT ATG AAA ATC GGC ATC ATT GGT GCA ATG G-3' (forward) and 5'-CTC **GCG** **GCC** **GCC** AGG CAA TCA CCA GAT CGG G-3' (reverse). Individual nucleotides that differed from the initial reported sequence of *pfs* are printed in bold. The nucleotide sequence has been submitted to the GenBank database (accession no. U24438).

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287

↳ coding sense strand

→ Double sequence → Reverse complement
→ down at stelle 5071

Met g s s H H s s g
5'-ATG GGC AGC AGC (CAT)₅ CAC AGC AGC ggc
L V P R G S H M
5'-CAT GAG ccg CGC GGC AGC CAT ATG - 3'

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E.coli MTA/SAH Nucleosidase 3.2.2.9

His-Tag Sequenz im Vektor (forward):

Thrombin recognition *NdeI*
 M G S₂ H₅ H S₂ G L V P | R G S H M
 5'-ATG GGC (AGC)₂ (CAT)₅ CAC (AGC)₂ GGC CTG GTG CCG CGC GGC AGC CAT ATG-3'

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Erkennungssequenz Thrombin → P₄-P₃-Pro-Arg/Lys-cut-P₁'-P₂' where P₃ and P₄ are hydrophobic and P₁' and P₂' are non-acidic

L V Q K L A H G *
 5'-CTG GTG CAG AAA CTT GCA CAT GGC TAA GTC ACT-3'
 3'-GAC CAC GTC TTT GAA CGT GTA CCG ATT
CTT AAG-5'
 Primer mit 5'-Überhang

Reverse primer:

EcoMTAN_EcoRI_rev

EcoRI
 5'-GAA TTC TTA GCC ATG TGC AAG TTT CTG CAC CAG-3'

Protease Schnittstellen:

Name of enzyme	No. of cleavages
Arg-C proteinase	6
Asp-N endopeptidase	12
Asp-N endopeptidase + N-terminal Glu	27
CNBr	4
Chymotrypsin-high specificity (C-term to [FYW], not before P)	11
Chymotrypsin-low specificity (C-term to [FYVML], not before P)	43
Clostrypain	6
Formic acid	12
Glutamyl endopeptidase	15
Hydroxylamine	2
LysC	11
LysN	11
NTCB (2-nitro-5-thiocyanobenzoic acid)	5
Pepsin (pH 1.3)	60
Pepsin (pH > 2)	53
Proline-endopeptidase [*]	1
Proteinase K	115
Staphylococcal peptidase I	13
Thermolysin	86
Trypsin	16

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Protocol I - CTAB protocol for the extraction of bacterial genomic DNA

This protocol is derived from the "miniprep" method described by Wilson (Wilson, K (1987) Preparation of genomic DNA from bacteria. In: Ausubel FM, Brent R, Kingston RE, Moore DD, Seidman JG, Smith JA, Struhl K, (eds) Current Protocols in Molecular Biology, pp. 2.4.1-2.4.5. John Wiley & Sons, New York). Broth cultures (2-5 ml) grown to mid-log growth phase are harvested in 2.0 ml Eppendorf tubes by centrifugation in a microfuge at 10,000-15,000 x g for 10-15 minutes. In general, late log growth phase cultures should not be used for preparing DNA, as nucleases tend to accumulate in older cultures. Alternatively, bacterial colonies grown on agar media may be washed off the agar and collected in an Eppendorf tube. The Bacteria Washing Buffer should not contain EDTA, as some bacteria (e.g., some Gram-negative species) will begin lysing upon exposure to chelating agents. After pelleting the cells, the medium is poured off and the rim of the tube is blotted with a paper towel to get rid of residual liquid. The bacterial pellet should weigh approximately 0.1 g (wet weight), which should provide 40-200 mg of DNA, depending upon the species of bacteria and the growth conditions. If there is more than 0.1 g per tube, the cell pellet should be resuspended with Bacteria Washing Buffer and redistributed accordingly into additional Eppendorf tubes. This is not unimportant, since the efficiency of the extraction decreases with increasing cell material. With experience, one can estimate reliably the mass of the cell pellet from its size.

0.5ml TN-Kultur $\xrightarrow{\text{Vor}}$ 2ml Zellpellet abfiltern
Steps in the protocol

1. Resuspend the cell pellet (approximately 0.1 g) completely with 564 μ l TE buffer (use a sterile ~~sterile~~ ^{sterile} pipette tip to mix the pellet and ensure complete resuspension).

(nur Kinnel)

2. Add approximately 10 μ g lysozyme (crystalline) to the cell suspension (from this point, do not vortex!). Mix thoroughly by inverting the Eppendorf tube several times. Incubate 10-60 minutes at 37°C. Add 6 μ l Proteinase K (10 mg/ml), and 30 μ l SDS (10-20%). Mix thoroughly (do not vortex!). Incubate at 37°C until the suspension becomes relatively clear and viscous.

3. Add 100 μ l NaCl (5 M) and mix thoroughly (do not vortex!). Incubate suspension at 65°C, 2 minutes. Add 80 μ l CTAB/NaCl solution (preheated at 65°C, use a pipette tip with the tip cut off to pipette the viscous CTAB/NaCl solution) and mix thoroughly (do not vortex!). Incubate suspension at 65°C, 10 minutes.

4. Extract suspension with an equal volume (approximately 800 μ l) chloroform:isoamyl alcohol (24:1) solution. Centrifuge (10,000 x g, 5 minutes). Transfer the upper (aqueous) phase (Supernatant 1), containing the nucleic acids, into a separate 2.0 ml Eppendorf tube.

5. Extract Supernatant 1 with an equal volume (approximately 800 μ l) of phenol/chloroform/isoamyl alcohol (25:24:1) solution. Centrifuge (15,000 x g, 5 minutes). Transfer the upper (aqueous) phase (Supernatant 2), containing the nucleic acids, into a separate 2.0 ml Eppendorf tube.

6. Extract Supernatant 2 with an equal volume (approximately 800 μ l) chloroform : isoamyl alcohol (24:1) solution. Centrifuge (10,000 x g, 5 minutes). Transfer the upper (aqueous) phase (Supernatant 3), containing the nucleic acids, into a separate 2.0 ml Eppendorf tube.

7. Add 0.7 volumes (approximately 560) isopropanol to precipitate nucleic acids. Mix gently by inverting the tube several times - the DNA should appear as a white, viscous, precipitate. Let sit at room temperature for 5 minutes to 1 hour. Centrifuge (12,000-15,000 x g, 15-30 minutes) at room temperature. The DNA should be visible as a pellet on the side of the Eppendorf tube. Remove the isopropanol carefully, so as to avoid disturbing the pellet.

8. Wash the pellet with 500 μ l EtOH (70%) by inverting the tube several times. Centrifuge 12,000-15,000 x g, 15-30 minutes at room temperature. Carefully remove the EtOH and blot the rim of the tube with a paper towel to get rid of excess liquid.

9. Briefly (not more than 5 minutes) dry pellet in a speed-vac. $5\text{ min ab in Kühlschrank}$

10. Resuspend each pellet in 50-60 μ l TE Buffer. Let sit at 37°C to allow the DNA to be resuspended completely.

11. Estimate the concentration of DNA in suspension by spectrophotometric measurement at 260 nm. For double-stranded DNA suspensions, at a wavelength of 260 nm and using a cuvette with a 1 cm light path, an OD of 1.0 is equal to a concentration of 50 mg/ml. The quality of the DNA can be estimated by measurements of the A260/A280 and the A260/A230 ratios. The size of the DNA can be estimated by agarose gel (0.5%, w/v) electrophoresis, subsequent staining with ethidium bromide and visualization by U.V. illumination. DNA of uniform size (approximately 20 kb) indicates that the DNA has been extracted without excessive shearing. DNA which has been sheared or degraded by nucleases will appear as a broad smear, or smaller molecular weight products.

$271.70 \text{ OD in TE-Buffer}$

Solutions

- Bacteria Washing Buffer 0.4 M NaCl, 50 mM Tris-HCl (pH 8.0), 50 mM MgSO₄, in sterile H₂O

- TE Buffer 10mMTris-HCl (pH 8.0), 1mM NaEDTA, in sterile H₂O

- NaCl 0.7 M NaCl in sterile H₂O

- CTAB/NaCl 10% (w/v) hexadecyltrimethyl ammonium bromide in sterile 0.7 M NaCl solution. [Heat solution to 65°C before bringing to final volume]

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+ 118,5 µl ddH₂O

2 Reaktionen


~~(1)~~ friner: EcoM~~X~~N-EcoR Γ _rev +
unel ~~(2)~~ EcoM~~X~~N-NdeI_FW

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ZY

- 10 g N-Z-mine AS (or my tryptic digest of casein, e.g. tryptone)
Sigma cat# N8317
- 5 g yeast extract
- We use the powdered extract from EM Sciences, available from TWR
TWR cat# E11 03733.0500
- 925 ml water

20x NPS

Component	100 ml	1 liter	seed liter
40 H ₂ O	90 ml	900 ml	-
(NH ₄) ₂ SO ₄	6.6 g	66 g	0.5 M
KH ₂ PO ₄	13.6 g	136 g	1 M
Na ₂ HPO ₄	14.2 g	142 g	1 M

- Add in sequence in beaker, stir until all dissolved
- pH of 20-fold dilution in water should be ~6.75

50x502: (5052 = 0.5 % glycerol, 0.05% glucose, 0.2% α -lactose)

Component	100 ml	1 liter
Glycerol (weigh in beaker)	25 g	250 g
H ₂ O	73 ml	730 ml
Glucose (Fisher FLBP254-1)	2.5 g	25 g
α -lactose (Sigma cat# L2425)	10 g	100 g

- Add in sequence in beaker, stir until all dissolved

- Lactose is slow to dissolve — may take two hours or more of stirring. Brief heating in a microwave can speed up dissolution of the lactose

1 M MgSO₄

- 24.65 g MgSO₄·7H₂O
- Water to make 100 ml

40% glucose (w/v)

Component	100 ml	200 ml
Glucose (Fisher FLBP254-1)	40 g	120 g
H ₂ O	74 ml	222 ml

- Add glucose to stirring water in beaker; DO NOT ATTEMPT TO ADD WATER TO GLUCOSE!
- Stir until all dissolved — may take 45 minutes or more of stirring.

80% glycerol (v/v) (= 100% w/v)

- 100 g glycerol (weigh in beaker)
- 20 ml water

20% α -lactose (w/v)

Component	100 ml	400 ml
α -lactose (Sigma cat# L2425)	20 g	120 g
H ₂ O	87.5 ml	525 ml

- add lactose to stirring water in beaker
- stir until all dissolved — may take 2 hours or more

1000x trace metals mixture (100 ml in ~50 mM HCl)

Prepare all metal stock solutions in dH₂O, except for FeCl₃, which is dissolved in ~0.1M HCl, as noted in the table below. Combine the metal solutions as in the table below.

Component	Vol	MW	1x conc
H ₂ O	36 ml	-	-
8.1 M FeCl ₃ ·6H ₂ O (dissolve in ~0.1 M HCl + 100 fold dil of cone HCl)	50 ml	270.30	50 μ M Fe
1M CuCl ₂	2 ml	110.99	20 μ M Cu
1M MnCl ₂ ·4H ₂ O	1 ml	197.91	10 μ M Mn
1M ZnCl ₂ ·H ₂ O	1 ml	287.56	10 μ M Zn
0.2 M CuCl ₂ ·2H ₂ O	1 ml	237.95	2 μ M Cu
0.1 M CuCl ₂ ·2H ₂ O	2 ml	170.486	2 μ M Cu
0.2 M NiCl ₂ ·6H ₂ O	1 ml	237.72	2 μ M Ni
0.1 M Na ₂ MoO ₄ ·2H ₂ O	2 ml	241.98	2 μ M Mo
0.1 M Na ₂ SeO ₃ ·5H ₂ O	2 ml	263.03	2 μ M Se
0.1 M H ₃ BO ₃	2 ml	61.83	2 μ M H ₃ BO ₃

Autoclave the stock solutions of the individual metals, except 0.1 M FeCl₃ in 1/100 volume cone HCl.A brief precipitate appeared upon addition of Na₂SeO₃, which redissolved rapidly.

Store at room temperature

When making growth media, add the metals mix before NPS. If NPS is already present when 1000x metals mix is added, a precipitate forms which disperses but retains a light turbidity. If the metals are diluted to near their final concentration before NPS is added, the medium remains clear. The metals also precipitate and disperse or redissolve when added to ZY; a precipitate caused by yeast extract. Although apparently not a problem, the precipitate could be avoided by diluting the metals in the water before dissolving the yeast extract in making ZY.

Antibiotics

- Kanamycin (25 mg/ml in water, filter sterilize)
- Chloramphenicol (25 mg/ml in 95% ethanol, filter sterilize)
- Ampicillin (50 mg/ml in water, filter sterilize)

P-0.5G defined medium for growth to saturation with little or no induction

- grow log-phase or saturated cultures for making feeder stocks and working stocks
- For all media, add 1 M MgSO₄ and 1000x metals mix before adding 20xNPS to avoid precipitate

Component	50 ml	100 ml	200 ml	Final concentration
40 H ₂ O (sterile)	-46.8 ml	-93.6 ml	-187.3 ml	-
1 M MgSO ₄	50 μ l	100 μ l	200 μ l	1 mM
10000x metals mix	5 μ l	10 μ l	20 μ l	1
40% glucose	0.625 ml	1.25 ml	2.5 ml	0.5%
20x NPS	2.5 ml	5 ml	10 ml	1 x
antibiotics	Choose only 1 from the list below			
kanamycin (25 mg/ml)	0.2 ml	0.4 ml	0.8 ml	100 μ g/ml
chloramphenicol (25 mg/ml)	50 μ l	100 μ l	200 μ l	25 μ g/ml

amp (50 mg/ml) 50 ml 100 ml 200 ml 50 mg/ml

ZYP-0.8G

- Rich medium for growth with little or no induction
- Cultures should go somewhat acid at saturation (slightly below pH 6)
- Collect cultures for freeze stocks well before saturation
- For all media, add 1 M MgSO₄ before adding 20nPS to avoid precipitate
- Kanamycin is used at significantly higher concentrations (100 µg/ml) than is normally the case (25-40 µg/ml). Shuler has found that in the T7 expression strains in these rich media, it does not provide adequate selection at the lower concentrations

Component	50 ml	100 ml	200 ml	400 ml	Concentration
ZY	-46.5 ml	-93 ml	-186 ml	-372 ml	-
1 M MgSO ₄	50 µl	100 µl	0.2 ml	0.4 ml	1 mM
40% glucose	1 ml	2 ml	4 ml	8 ml	0.8%
20nPS	2.5 ml	5 ml	10 ml	20 ml	1x
Antibiotics, as needed					
Choose only 1 from the list below					
kanamycin (25 µg/ml)	200 µl	0.4 ml	0.8 ml	1.6 ml	100 µg/ml
chloramphenicol (25 µg/ml)	50 µl	100 µl	0.2 ml	0.4 ml	25 µg/ml
ampicillin (50 µg/ml)	50 µl	100 µl	0.2 ml	0.4 ml	50 µg/ml

+ 0.2x
Trace
Metals

ZYP-5052 rich medium for auto-induction

- For all media, add 1 M MgSO₄ before adding 20nPS to avoid precipitate
- Kanamycin is used at significantly higher concentrations (100 µg/ml) than is normally the case (25-40 µg/ml). Shuler has found that in the T7 expression strains in these rich media, it does not provide adequate selection at the lower concentrations
- We use 400 ml in a 2 liter baffled flask

Adequate aeration is essential to the performance of this media. Don't even think about using more than 20% of the nominal volume of the flask. Baffled flasks will give significantly better performance. You might obtain adequate results with non-baffled flasks, but I don't recommend it.

Component	200 ml	400 ml	500 ml	1 liter	Concentration
ZY	-186 ml	-372 ml	-464 ml	-928 ml	-
1 M MgSO ₄	0.2 ml	0.4 ml	0.5 ml	1 ml	1 mM
Sta's 5452	4 ml	8 ml	10 ml	20 ml	1x
20nPS	10 ml	20 ml	25 ml	50 ml	1x
Antibiotics					
Choose only 1 from the list below					
kanamycin (25 µg/ml)	0.8 ml	1.6 ml	2 ml	4 ml	100 µg/ml
chloramphenicol (25 µg/ml)	0.2 ml	0.4 ml	0.5 ml	1 ml	25 µg/ml
ampicillin (50 µg/ml)	0.2 ml	0.2 ml	0.5 ml	1 ml	50 µg/ml

+ 0.2x
Trace
Metals

LB Agar plates + 1% glucose

- Sterile LB agar as usual for pouring plates. Allow agar to cool to 55-60 °C, usually in a 56 °C water bath.
 - Add 1/40 volume of 40% glucose (e.g. for 400 ml, add 1 ml glucose)
 - Add antibiotic to an appropriate concentration
 - Pour the plate

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
30.05.2012 16:33:31

Enhancement: 1
Exposure Time: 100Brightness: 50
Contrast: 50

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Lysis of Cells

9. Allow the frozen bacterial cell pellet from Step 7 to thaw at room temperature for 5–10 minutes. Resuspend the pellet in 18 ml (10 ml) of Alkaline lysis solution I.

The volumes given in parentheses in the remainder of this protocol should be used only with cultures that have been treated with chloramphenicol at Step 4.

10. Add 2 ml (1 ml) of a freshly prepared solution of 10 mg/ml lysozyme.

11. Add 40 ml (20 ml) of freshly prepared Alkaline lysis solution II. Close the top of the centrifuge bottle and mix the contents thoroughly by gently inverting the bottle several times. Incubate the bottle for 5–10 minutes at room temperature.

Protocol note: Some of the superheated alkali results in reversible denaturation (Vinegar and Lebowitz 1986). The resulting cyclic coiled DNA cannot be cleaved with restriction enzymes, and it migrates through agarose gels about twice the rate of superhelical DNA and stains poorly with ethidium bromide. Traces of this form of DNA can be seen in plasmids prepared by alkaline lysis of bacteria.

12. Add 20 ml (15 ml) of ice-cold Alkaline lysis solution III. Close the top of the centrifuge bottle and mix the contents gently but well by swirling the bottle several times (there should no longer be two distinguishable liquid phases). Place the bottle on ice for 10 minutes.

A flocculent white precipitate consisting of chromosomal DNA, high-molecular-weight RNA, and potassium (S1)3/protein/cell wall complexes will form during this incubation.

Potassium acetate is used in preference to sodium acetate in Alkaline lysis solution III because the potassium salt of dodecyl sulfate is far less soluble than the sodium salt.

13. Centrifuge the bacterial lysate at ≥20,000g (11,000 rpm in a Sorvall GSA rotor) for 30 minutes at 4°C in a medium-speed centrifuge. Allow the rotor to stop without braking. At the end of the centrifugation step, decant the clear supernatant into a graduated cylinder. Discard the pellet remaining in the centrifuge bottle.

The failure to form a compact pellet after centrifugation is usually a consequence of inadequate mixing of the bacterial lysate with Alkaline lysis solution II (Step 11). If the bacterial debris does not form a tightly packed pellet, centrifuge again at 20,000g (11,000 rpm in a Sorvall GSA rotor) for a further 15 minutes, and then transfer as much of the supernatant as possible to a fresh tube. Decanting the supernatant through four-ply gauze at this step helps to remove traces of viscous genomic DNA and protein precipitate.

Recovery of Plasmid DNA

14. Measure the volume of the supernatant. Transfer the supernatant together with 0.6 volume of isopropanol to a fresh centrifuge bottle. Mix the contents well and store the bottle for 10 minutes at room temperature.

15. Recover the precipitated nucleic acids by centrifugation at 12,000g (8000 rpm in a Sorvall GSA rotor) for 15 minutes at room temperature.

Salt may precipitate if centrifugation is carried out at 4°C.

16. Decant the supernatant carefully, and invert the open bottle on a paper towel to allow the last drops of supernatant to drain away. Rinse the pellet and the walls of the bottle with 70% ethanol at room temperature. Drain off the ethanol, and use a Pasteur pipette attached to a

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dilute
 Pellets

Wesentlich
 kompaktere
 Pellets

Protocol 3: Preparation of Plasmid DNA by Alkaline Lysis with SDS: Maxipreparation 1.41

vacuum line to remove any beads of liquid that adhere to the walls of the bottle. Place the inverted, open bottle on a pad of paper towels for a few minutes at room temperature to allow the ethanol to evaporate.

If the pellet of DNA is dried in a desiccator or under vacuum, it becomes difficult to dissolve under some circumstances and may denature (Swane et al. 1987). Drying the pellet for 10–15 minutes at room temperature is usually sufficient for the ethanol to evaporate without DNA becoming dehydrated.

17. Dissolve the damp pellet of nucleic acid in 3 ml of TE (pH 8.0).

18. Purify the crude plasmid DNA either by column chromatography (Protocol 9), precipitation with polyethylene glycol (Protocol 8), or equilibrium centrifugation in CsCl-ethidium bromide gradients (Protocols 10 and 11).

19. Check the structure of the plasmid by restriction enzyme digestion followed by gel electrophoresis.

For recommendations on troubleshooting, please see Table 1-5.

The typical yield of high-copy-number plasmid vectors or of amplified low-copy-number vectors prepared by this method is ~3–5 µg of DNA/g of original bacterial culture. The yield of recombinant plasmid DNA is usually lower than the yield of the original vector because of the loss of the cloned DNA fragment. Yields of <10 µg/ml indicate that the plasmid DNA is toxic to *E. coli*; the complete absence of DNA indicates that the plasmid DNA is toxic to *E. coli* during extraction and purification. The first problem can often be solved by switching to a low-copy-number vector or to a vector that carries prokaryotic transcriptional termination signals. The second problem is often the result of accidental loss of the plasmid DNA pellet after ethanol precipitation. It also occurs when the recombinant plasmid is large enough to be precipitated together with the chromosomal DNA after addition of potassium acetate. Such a large plasmid will also be susceptible to breakage during extraction and purification. This difficulty can be solved by using a more gentle procedure of cell lysis, such as that described in Protocol 7.

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METHOD

- Transfer 3 ml of the crude large-scale plasmid preparation to a 15-ml Corex tube and chill the solution to 0°C in an ice bath.
- Add 3 ml of an ice-cold solution of 5 M LiCl to the crude plasmid preparation, mix well, and centrifuge the solution at 12,000g (10,000 rpm in a Sorvall SS-34 rotor) for 10 minutes at 4°C.
- Transfer the supernatant to a fresh 30-ml Corex tube. Add an equal volume of isopropanol. Mix well. Recover the precipitated nucleic acids by centrifugation at 12,000g (10,000 rpm in a Sorvall SS-34 rotor) for 10 minutes at room temperature.
- Decant the supernatant carefully, and invert the open tube to allow the last drops of supernatant to drain away. Rinse the pellet and the walls of the tube with 70% ethanol at room temperature. Carefully discard the bulk of the ethanol, and then use a vacuum aspirator to remove any beads of liquid that adhere to the walls of the tube. Place the inverted, open tube on a pad of paper towels for a few minutes until no trace of ethanol is visible. At this stage, the pellet should still be damp.
- Dissolve the damp pellet of nucleic acid in 500 µl of TE (pH 8.0) containing RNase A. Transfer the solution to a microfuge tube and store it for 30 minutes at room temperature.

Protocol 8: Purification of Plasmid DNA by Precipitation with Polyethylene Glycol 1.61

- Extract the plasmid-RNase mixture once with phenol:chloroform and once with chloroform.
- Recover the DNA by standard ethanol precipitation.
- Dissolve the pellet of plasmid DNA in 1 ml of sterile H₂O, and then add 0.5 ml of PEG-MgCl₂.
- Store the solution for ≥10 minutes at room temperature, and then collect the precipitated plasmid DNA by centrifugation at maximum speed for 20 minutes at room temperature in a microfuge.
- Remove traces of PEG by resuspending the pellet of nucleic acid in 0.5 ml of 70% ethanol. Collect the nucleic acid by centrifugation at maximum speed for 5 minutes in a microfuge.
- Remove the ethanol by aspiration and repeat Step 10. Following the second rinse, store the open tube on the bench for 10–20 minutes to allow the ethanol to evaporate.
- Dissolve the damp pellet in 500 µl of TE (pH 8.0). Measure the OD₂₆₀ of a 1:100 dilution in TE (pH 8.0) of the solution, and calculate the concentration of the plasmid DNA assuming that 1 OD₂₆₀ = 50 µg of plasmid DNA/ml.
- For information on absorption spectroscopy of DNA, please see Appendix 8.
- Store the DNA in aliquots at ~20°C.

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
05.06.2012 12:01:40

Enhancement: 110
Exposure Time: 323Brightness: 50
Contrast: 50
Gamma: 25

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
S-Gal - Agar

S-Gal™

Product Code S 9811
Store at Room Temperature

Product Information

Precautions and Disclaimer

Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

When dry-blended with isopropyl β -D-thiogalactoside in agar medium, S-Gal is autoclavable. There's no waiting time for media to cool before adding the substrate. Additionally, since it can be blended with growth medium prior to autoclaving, making up stock solutions in dimethylformamide (DMF) or dimethyl sulfoxide (DMSO) is not necessary. Should you desire to work from a concentrated solution, S-Gal is soluble in DMSO at 50 mg/ml.

S-Gal is added to the medium prior to autoclaving at a recommended concentration of 300 mg/l of medium, along with 500 mg/l of ferrous ammonium citrate.

Note: The iron(III) (feric or Fe^{3+}) ion is required for color development and must be added to any S-Gal formulation. Pre-formulated products such as C 4478, S-Gal™/LB Agar Blend contain ferrous ammonium citrate (F 5879). A medium prepared with S-Gal is moderately dark due to the presence of ferrous ammonium citrate. This darker background often provides enhanced contrast for automated colony counting or isolation.

Storage/Stability

The product as supplied is stored at room temperature. In prepared medium, it is stable for two weeks when stored at 4 °C.

References

- U.S. Patent #6,008,008.
- Heuermann, K. and Cosgrove, J., S-Gal™: A superior dye to X-gal for clonal selection. LifeScience Quarterly, 2(2), 2-4 (2001).
- Heuermann, K. and Cosgrove, J., An autoclavable dye for color selection of cloned DNA inserts. BioTechniques, 30(5), 1142-1147 (2001).

GWS/MAM 11/01

Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

2x 100 ml L-agar:

10 g/l Trypton / (Pepton)

5 g/l Hefeextrakt

5 g/l NaCl



1,5% (w/v) Agar

500 mg/l Ammonium eisen(II)-citrat

300 mg/l S-gal

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Bacto - XTT - Agar

- last experiment in molecular genetics

Difco Antibiotic medium #2:

Beef Extract
Yeast Extract
Pepton
Agar

15g/l
3g/l
6g/l
15g/l

} 25,5g medium + 950ml
ddH₂O (TCC)
- 50mg 2,3,5-triphenyltetrazolium
- 50ml 20% Bactose

→ autoclave

Eckmechan:

- 1) 1,5g Bacto beef extract, 3g yeast extract, 6g Pepton, 15g agar
- 2) 23g Difco nutrient agar, 8g NaCl
pro Liter

Änderung

- 6g Pepton Trypton

- 3g Hefeextrakt

- 8,5g Agar

- 50ml 20% (w/v) Bactose


- 50mg XTT (2,3,5-triphenyl-tetrazolium-5-sulfophenyl)

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After 24 hours at 37°C. a - *E. coli* DH5 α with pMLB1034 vector. b - *E. coli* BL21(DE3) with pMLB1034 Isr'-LacZ $^+$ transcriptional fusion. c - *E. coli* DH5 α with pMLB1034 Isr'-LacZ $^+$ transcriptional fusion.

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Protocol 3

Preparation of Plasmid DNA by Alkaline Lysis with SDS: Maxipreparation

PLASMID DNA MAY BE ISOLATED FROM LARGE-SCALE (500 ml) bacterial cultures by treatment with alkali and SDS. The resulting DNA preparation may be further purified by column chromatography or centrifugation through CsCl-ethidium bromide gradients.

MATERIALS

CAUTION: Please see Appendix 12 for appropriate handling of materials marked with <1>.

Buffers and Solutions

Please see Appendix 1 for components of stock solutions, buffers, and reagents.

Alkaline lysis solution I

For preparations of plasmid DNA that are to be subjected to further purification by chromatography (please see Protocol 9), sterile Alkaline lysis solution I may be supplemented just before use with the appropriate volume of 20 mg/ml case-free RNase A (pancreatic RNase) to give a final concentration of 100 µg/ml. Addition of RNase is not recommended at this stage if the DNA is to be further purified by other methods (please see Protocols 10 and 11).

Alkaline lysis solution II

Solution II should be freshly prepared and used at room temperature.

Alkaline lysis solution III

Antibiotic for plasmid selection

Chloramphenicol (34 mg/ml) <1>

Optional, please see Step 4.

Ethanol

Isopropanol

STE

TE (pH 8.0)

Enzymes and Buffers

Lysozyme (10 mg/ml)

Please see the information panel on LYSOZYMES.

Restriction endonucleases

Gels

Agarose gels

Please see Steps 8 and 19.

Protocol 3: Preparation of Plasmid DNA by Alkaline Lysis with SDS: Maxipreparation 1.39

Media

LB, YT, or Terrific Broth

Centrifuges and Rotors

Sorvall CSA rotor or equivalent

Sorvall SS-34 rotor or equivalent

Additional Reagents

Steps 8 and 19 of this protocol require reagents listed in Chapter 5, Protocol 1.

Step 18 of this protocol requires reagents listed in Protocol 8, 9, 10, or 11 of this chapter.

METHOD

Preparation of Cell Culture

- Inoculate 30 ml of rich medium (LB, YT, or Terrific Broth) containing the appropriate antibiotic either with a single colony of transformed bacteria or with 0.1–1.0 ml of a small-scale liquid culture grown from a single colony.

To ensure that the culture is adequately aerated:

- The volume of the culture flask should be at least four times greater than the volume of the bacterial culture.
- The culture flask should be loosely capped.
- The culture should be incubated with vigorous agitation.

- Incubate the culture at the appropriate temperature with vigorous shaking until the bacteria reach late log phase ($OD_{600} = -0.6$).

- Inoculate 500 ml of LB, YT, or Terrific Broth medium (prewarmed to 37°C) containing the appropriate antibiotic in a 2-liter flask with 25 ml of the late-log-phase culture. Incubate the culture for ~2.5 hours at 37°C with vigorous shaking (300 cycles/minute on a rotary shaker).

The OD_{600} of the resulting culture should be ~0.4. Because the growth rates of different bacterial strains will vary, the culture may have to be incubated slightly longer or shorter than 2.5 hours to reach an OD of 0.4.

- For relaxed plasmids with low or moderate copy numbers, add 2.5 ml of 34 mg/ml chloramphenicol solution. The final concentration of chloramphenicol in the culture should be 170 µg/ml.

▲ IMPORTANT For high-copy-number plasmids, do not add chloramphenicol.

- Incubate the culture for a further 12–16 hours at 37°C with vigorous shaking (300 cycles/minute on a rotary shaker).

- Remove an aliquot (1–2 ml) of the bacterial culture to a fresh microfuge tube and store it at 4°C. Harvest the remainder of the bacterial cells from the 500-ml culture by centrifugation at 2700g (4100 rpm in a Sorvall GSA rotor) for 15 minutes at 4°C. Discard the supernatant. Stand the open centrifuge bottle in an inverted position to allow all of the supernatant to drain away.

- Resuspend the bacterial pellet in 200 ml of ice-cold STE. Collect the bacterial cells by centrifugation as described in Step 6. Store the pellet of bacteria in the centrifuge bottle at ~20°C.

- Use one of the methods described in Protocol 1 or Protocol 4 to prepare plasmid DNA from the 1–2-ml aliquot of bacterial culture set aside in Step 6. Analyze the minipreparation plas-

Lysis of Cells

mid DNA by digestion with the appropriate restriction enzyme(s) and agarose gel electrophoresis to ensure that the correct plasmid has been propagated in the large-scale culture. This kind of control may seem a little overly compulsive. However, it provides valuable insurance against errors that may be difficult to retrieve and may cause considerable loss of time.

9. Allow the frozen bacterial cell pellet from Step 7 to thaw at room temperature for 5–10 minutes. Resuspend the pellet in 18 ml (10 ml) of Alkaline lysis solution I. The volumes given in parentheses in the remainder of this protocol should be used only with cultures that have been treated with chloramphenicol at Step 4.
10. Add 2 ml (1 ml) of a freshly prepared solution of 10 mg/ml lysozyme.
11. Add 40 ml (20 ml) of freshly prepared Alkaline lysis solution II. Close the top of the centrifuge bottle and mix the contents thoroughly by gently inverting the bottle several times. Incubate the bottle for 5–10 minutes at room temperature. Prolonged exposure of superhelical DNA to alkali results in irreversible denaturation (Vingrad and Linn 1966). The resulting cyclic coiled DNA cannot be cleaved with restriction enzymes, and it migrates through agarose gels at about twice the rate of superhelical DNA and stains poorly with ethidium bromide. Traces of this form of DNA can be seen in plasmids prepared by alkaline lysis of bacteria.
12. Add 20 ml (15 ml) of ice-cold Alkaline lysis solution III. Close the top of the centrifuge bottle and mix the contents gently but well by swirling the bottle several times (there should be two distinguishable liquid phases). Place the bottle on ice for 10 minutes. A flocculent white precipitate consisting of chromosomal DNA, high-molecular-weight RNA, and potassium/S1/protein/cell wall complexes will form during this incubation. Potassium acetate is used in preference to sodium acetate in Alkaline lysis solution III because the potassium salt of dodecyl sulfate is far less soluble than the sodium salt.
13. Centrifuge the bacterial lysate at $\geq 20,000g$ (11,000 rpm in a Sorvall GSA rotor) for 30 minutes at 4°C in a medium-speed centrifuge. Allow the rotor to stop without braking. At the end of the centrifugation step, decant the clear supernatant into a graduated cylinder. Discard the pellet remaining in the centrifuge bottle. The failure to form a compact pellet after centrifugation is usually a consequence of inadequate mixing of the bacterial lysate with Alkaline lysis solution II (Step 11). If the bacterial debris does not form a tightly packed pellet, centrifuge again at 20,000g (11,000 rpm in a Sorvall GSA rotor) for a further 15 minutes, and then transfer as much of the supernatant as possible to a fresh tube. Decanting the supernatant through four-ply gauze at this step helps to remove traces of viscous genomic DNA and protein precipitate.

Recovery of Plasmid DNA

14. Measure the volume of the supernatant. Transfer the supernatant together with 0.6 volume of isopropanol to a fresh centrifuge bottle. Mix the contents well and store the bottle for 10 minutes at room temperature.
15. Recover the precipitated nucleic acids by centrifugation at 12,000g (8000 rpm in a Sorvall GSA rotor) for 15 minutes at room temperature. Salt may precipitate if centrifugation is carried out at 4°C.
16. Decant the supernatant carefully, and invert the open bottle on a paper towel to allow the last drops of supernatant to drain away. Rinse the pellet and the walls of the bottle with 70% ethanol at room temperature. Drain off the ethanol, and use a Pasteur pipette attached to a

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17. Dissolve the damp pellet of nucleic acid in 3 ml of TE (pH 8.0).
18. Purify the crude plasmid DNA either by column chromatography (Protocol 9), precipitation with polyethylene glycol (Protocol 8), or equilibrium centrifugation in CsCl-ethidium bromide gradients (Protocols 10 and 11).
19. Check the structure of the plasmid by restriction enzyme digestion followed by gel electrophoresis.

For recommendations or troubleshooting, please see Table 1-5.

The typical yield of high-copy-number plasmid vectors or of amplified low-copy-number vectors prepared by alkaline lysis is 3–5 µg of DNA/ml of original bacterial culture. The yield of recombinant plasmids containing inserted foreign DNA is usually slightly lower, depending on the size and nature of the cloned DNA fragment. Yields of < 10 µg/ml indicate that the plasmid DNA is toxic to *E. coli*; the complete absence of DNA indicates that the plasmid has been lost during extraction and purification. The first problem can often be solved by switching to a lower-copy-number vector and/or to a vector that carries prokaryotic transcriptional termination signals. The second problem is often the result of accidental loss of the plasmid DNA pellet after ethanol precipitation. It also occurs when the recombinant plasmid is large enough to be precipitated together with the chromosomal DNA after addition of potassium acetate. Such a large plasmid will also be susceptible to breakage during extraction and purification. This difficulty can be solved by using a more gentle procedure of cell lysis, such as that described in Protocol 7.

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