

WEB308 - SOMT ammonium sulfate precipitation

Benjamin Weigel

23rd September, 2014

Introduction

- refolded 2.5 mL 1 mg/mL SOMT in solubilization buffer in 50 mL Buffer 12 by fast dilution (in 50 mL tube)
- incubation over night at 4°C
- centrifuged at 16.000 $x g$ / 4°C to remove precipitate
- took 5 mL aliquots and performed ammonium sulfate precipitation test (see Figure 1)




Figure 1: Experimental setup for the ammonium sulfate precipitation.

Table 1: Amount of Ammonium sulfate to add to 5 mL solution at each step.

AS cut	20-30%	30-40%	40-50%	50-60%	60-70%
AS (1 st step)	0.53 g	0.82 g	1.13 g	1.455 g	1.8 g
AS (2 nd step)	0.275 g	0.28 g	0.29 g	0.3 g	0.31 g

1 SDS-PAGE

- ca. 50 $\mu\text{g/mL}$ protein in refolding reaction
- precipitation via ammonium sulfate
- SDS-PAGE analysis of pellet & supernatant of first step
- SDS-PAGE analysis of pellet & supernatant of second step \rightarrow 16 samples
- 200 μL = 10 μg protein \rightarrow take 200 μL from precipitation in eppi & centrifuge
- resuspend AS pellet (1) in 50 μL PBS \rightarrow 5 μL for Gel
- TCA precipitation with supernatant, take up pellet from TCA-precipitation in 10 μL PBS and analyze via SDS-PAGE
- resuspend protein pellet from AS-precipitation (2) in 0.25 mL PBS \rightarrow 10 μL (ca 10 μg) on gel




Figure 2: Ammonium sulfate precipitation test straight from refolding SOMT in refolding buffer 12. The majority of SOMT protein is precipitated in the 30–40% cut of ammonium sulfate saturation (Fig. 2b). However, there is a some SOMT still in solution even at a 50% AS saturation (Fig. 2a).

Table 2: Setup of the gels for SDS-PAGE analysis

lane	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
pellets	M	20%	30%	40%	50%	M	30%	40%	50%	60%	M					
			1 st precipitation					2 nd precipitation								
supernatant	M	20%	30%	40%	50%	M	30%	40%	50%	60%	M	MTAN				
			1 st precipitation					2 nd precipitation								